



	<p>UVEX Respirator mask FFP3</p>	
	<p>Area of application and properties</p>	
<ul style="list-style-type: none"> • Maximum protection and comfortable to wear: Half-mask with a high filtering capacity (> 99%). The innovative design and optimised shape of this FFP3 preformed mask together with the soft sealing lip means it is extremely comfortable to wear, as well as offering optimized protective functions. The soft material and the special shape allow for a perfect fit. Minimal respiratory resistance and the exhalation valve also increase the comfort factor. Type tested and certified in accordance with EC directive 89/686/EEC (Personal protective equipment). • Application area: Protection from solid and liquid particles, e.g. for handling CMR¹⁾ drugs (e.g. cytostatics, virostatics, biological agents and viruses). • Protection capacity: Protection from all CMR drugs or chemicals is not guaranteed. Does not protect against gases and vapours. • Replacement interval: Depends on inhalation resistance. Maximum use for one shift (8 hours) with at least three breaks in between. For use with CMR-drugs, a maximum of 4 hrs, as long as there is no obvious contamination. Replace immediately in the event of visible contamination. Single use – Do not re-use respirator mask! • Before use: Check for any damage. Do not use damaged masks. • Limitation: Single use mask! Laboratory test results can differ from values obtained in practice considerably and, as a result, the useful life of the mask may increase or decrease. Facial hair and scars can impair the fit of the mask and therefore the level of protection. • Disposal: If used in combination with CMR-drugs and contaminated areas/items: waste requiring supervision (waste code: 18 01 04 in accordance with 2000/532/EC); in the event of heavy contamination, waste requiring special supervision (waste code: 18 01 08* in accordance with 2000/532/EC); collect and dispose of waste separately! • ¹⁾: Carcinogenic; mutagenic; toxic to reproduction 		
<p>Type</p>		
Size	Universal	
Item no.	121287	
Packaging unit (PU)	1	
Dimensions of single PU	140 mm x 160 mm	
Colour of mask	White	
Colour of valve	Silver valve	
<p>Material properties</p>		
Headband	Thermoplastic elastomer	
Weight	Excluding packaging: 25 g	
<p>Particle filtration efficiency</p>		
Testaerosol	Minimum filtration efficiency	
Sodium chloride	> 99%	
Paraffin oil	> 99%	

Dolomite clogging test (for FFP3-mask) (Dolomite DRB4/15b test aerosol)	Passed
Gases	Unsuitable
CE-mark	
CE-mark according to the PPE directive 89/686/EEC for complex PPE of the category III . The type-tested certificate is based on DIN EN 149:2001. Documented by the EC type test certificate of the manufacturer. Quality assurance (EC-quality assurance system with monitoring): control measures in compliance with 89/686/EEC by the applicable notified body CIOP PIB (1437).	
Notified body „1437“	
CIOP PIB (Zentralinstitut für Arbeitsschutz – Nationales Forschungsinstitut) Prüf-, Überwachungs- und Zertifizierungsstelle Nr. 1437 , Czerniakowska 16, 00-701 Warsaw, Poland.	
	Quality management system
Our quality management system is tested and certified by the TÜV Product Service GmbH (a certification body by the German Accreditation Council) in accordance with DIN EN ISO 9001:2008. Regular audits and production site inspections guarantee the quality of our products.	
	Storage and transport conditions
<ul style="list-style-type: none"> • Dark (protect from direct UV and sunlight) • Cool (-20 to +40°C) • Dry • No contact with pointed and/or sharp objects 	
	Shelf life
<ul style="list-style-type: none"> • 3 years 	
Manufacturer	
UVEX Arbeitsschutz GmbH, Würzburger Str. 181 – 189, 90766 Fürth, Germany	
Distributor	
BERNER INTERNATIONAL GmbH Mühlenkamp 6 - 25337 Elmshorn - Germany Tel.: +49/(0)4121 43560 - Fax: +49/(0)4121 4356 20 Info@berner-international.de www.berner-international.eu	