

# BIOCLEAN PRELUDE™ Sterile Latex Gloves BPSL

## About

Designed for easy double-donning the BioClean Prelude™ Sterile Latex Cleanroom Gloves also provide exceptional flexibility and wearer comfort. The natural-coloured 300mm (12") long, latex gloves feature a textured surface for enhanced grip and a beaded cuff for stability on the arm.

## Specifications

**COMPATIBILITY:**

ISO Class 4 & EU GMP Grade A

**LENGTH:**

300mm (12")

**MATERIAL:**

Latex

**PROTEIN LEVEL:**

Less than 50 µg/gram

**SURFACE:**

Textured

**SHAPE:**

Ambidextrous

**COLOUR:**

Natural

**STERILIZATION:**

Gamma irradiation


## Features

- Powder-free
- Beaded cuff
- EasyOn™ proprietary technology allows for easy double-donning
- Non-particulating EasyTear™ packaging

TEXTURED FOR IMPROVED GRIP

EXCEPTIONAL FLEXIBILITY & COMFORT

DOUBLE-DONNABLE



Ms. G. Love™ represents BioClean's range of cleanroom gloves

## Quality Standards

- Manufactured in a facility holding ISO 9001: 2008 and 89/686/EEC Article 11b certifications
- Conforms to Category 3 - Complex Design Personal Protective Equipment PPE Directive 89/686/EEC
- Processed in an NEBB certified ISO Class 4 environment
- Complies with the requirements of Glove Standard EN420:2003 and EN374-1, 2 & 3:2003
- Physical properties comply with European Medical Glove Standard EN455-2:2009
- Sterilized in accordance with ISO 11137-1:2006

## Physical Properties

CHARACTERISTICS	VALUE	TEST METHOD
Freedom From Holes	1.5 AQL Performance Level 2	EN 374-2:2003
Minimum Length (mm/inches)	290/11.4	EN 420:2003
Single Wall Thickness Minimum (mm)	Cuff: 0.11 Palm: 0.17 Finger Tip: 0.20	EN 455-2:2009
Force At Break During Shelf Life	≥ 9 N	EN 455-2:2009

PALM WIDTHS (mm) ± 5mm						
SIZE	5.0-5.5	6.0-6.5	7.0-7.5	8.0-8.5	9.0	10.0
mm	77	83	95	106	112	118
TEST METHOD				EN455-2:2009		

## Sterilization

Method: Gamma irradiation  
Minimum Dose: 25kGy  
Sterility Assurance Level: 10<sup>-6</sup>

## Shelf Life & Storage

Five (5) years from date of manufacture.  
Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## Glove Cleanliness Characteristics

PARTICLES				TYPICAL PARTICLE COUNT				TEST METHOD
≥ 0.5µm(counts/cm <sup>2</sup> )				<1500				IEST-RP-CC005.3
TYPICAL EXTRACTABLE IONS								TEST METHOD
ANIONS	Fluoride	Chloride	Nitrite	Bromide	Nitrate	Phosphate	Sulphate	IEST-RP-CC005.3
(µg/cm <sup>2</sup> )	ND	0.97	0.67	NT	0.67	0.07	ND	
CATIONS	Lithium	Sodium	Ammonium	Potassium	Calcium	Magnesium	Zinc	
(µg/cm <sup>2</sup> )	NT	0.02	NT	0.03	0.37	0.003	ND	
ND = Not Detected, NT = Not Tested								

## TO ORDER

Re-order code	Size	Packing
BPSL-5055	5.0-5.5	One pair per inner PE wallet; one wallet per sealed EasyTear™ PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs)
BPSL-6065	6.0-6.5	
BPSL-7075	7.0-7.5	
BPSL-8085	8.0-8.5	
BPSL-9090	9.0	
BPSL-1010	10.0	

## CONTACT

### Europe (Headquarters)

Nitritex Ltd.  
Minton Enterprise Park,  
Oaks Drive,  
Newmarket, Suffolk,  
CB8 7YY,  
United Kingdom  
T: +44 (0) 1638 663338  
F: +44 (0) 1638 668890  
E: [europa@bioclean.com](mailto:europa@bioclean.com)

### Other Regions

Asia: [malaysia@bioclean.com](mailto:malaysia@bioclean.com)  
Americas: [americas@bioclean.com](mailto:americas@bioclean.com)



BioClean™ is a brand of Nitritex Ltd.  
For more information, please visit [www.bioclean.com](http://www.bioclean.com)

BPSL/1602/PDS4