

SPORICIDE
LOW
ENDOTOXINS
 NO RESIDUE
 GOLD STANDARD
 6% HYDROGEN PEROXIDE
 NO RESIDUE
 TRIPLE BAGGED
 SPORICIDE
 NO RESIDUE
IRRADIATED
TRIPLE BAGGED
 SPORICIDE
 TRIPLE BAGGED
 NO RESIDUE
IRRADIATED
 TERMINAL STERILISATION
DISINFECTION
WATER FOR INJECTION
 LOW ENDOTOXINS
 GOLD STANDARD
 TERMINAL 6% HYDROGEN PEROXIDE
 STERILISATION
DISINFECTION
WATER FOR INJECTION
IRRADIATED
GOLD STANDARD
SPORICIDE
 NO RESIDUE
 6% HYDROGEN PEROXIDE
 PRODUCT CONTACT AREAS
GOLD STANDARD
IRRADIATED
 6% HYDROGEN PEROXIDE
 LOW ENDOTOXINS
NO RESIDUE
IRRADIATED
 6% HYDROGEN PEROXIDE
 LOW ENDOTOXINS
NO RESIDUE
TRIPLE BAGGED
WATER FOR INJECTION
HYDROPURE
 DISINFECTION
 SPORICIDE



6%
HYDROGEN PEROXIDE
SPORICIDAL
 WITH NO RESIDUE


CONTEC
www.contecinc.com

Contec HydroPure is a sporocidal biocide which leaves no residue and has a guaranteed endotoxin limit of <0.25 EU/ml. 6% stabilised peroxide blended with Water for Injection **Contec Sterile HydroPure** is specifically designed for product contact areas. The only terminally sterilised hydrogen peroxide available.

CONTEC HYDROPURE

THE WORLD'S FIRST IRRADIATED H₂O₂ PRODUCT

Contec HydroPure is specifically designed for product contact areas in high grade cleanrooms. Manufactured in a Grade C cleanroom and filtered to 0.2 µm, Contec HydroPure has no residue and a guaranteed endotoxin limit of ≤0.25EU/ml. Contec HydroPure is provided sterile, via gamma irradiation at no less than 25 kGy.

WHAT IS IT ?

Contec HydroPure is a sterile ready to use sporicidal biocide containing 6% v/v hydrogen peroxide in water for injection. It has a guaranteed endotoxin level of less than 0.25EU/ml and is triple bagged for ease of transfer into high grade cleanrooms.

WHY IS IT DIFFERENT ?

Contec has been able to create a hydrogen peroxide which is suitable for terminal sterilisation. All other commercially available sterile hydrogen peroxide products are rendered sterile by either a sterile or aseptic fill process. GMP suggests that wherever terminal sterilisation in the final container is possible it should be employed.

WHY IS IT SO SUITABLE ?

The control of residues and endotoxins is key in product contact areas, Contec HydroPure breaks down into water and oxygen leaving no residue behind on the surface. Since the presence of bacterial spores in these areas is negligible, there is less requirement for a very powerful sporicide, Contec HydroPure gives a log 3 reduction in spores in 15 mins *. Contec HydroPure is blended with water for injection and then filtered to 0.2µm in a Grade C cleanroom to minimise any bioburden in the product allowing a product to be created with a guaranteed endotoxin limit of ≤0.25EU/ml.



WHAT DO YOU GET ?

A terminally sterilized, ready to use, disinfectant which is :-

Sporicidal – sporicidal in 15 mins *

Endotoxin controlled – less than 0.25EU/ml

Residue free – breaks down to water and oxygen on a surface

Cleanroom compatible – filtered and triple bagged in a cleanroom environment

Fully validated – provided with complete testing package to EN norms

For further information or to request a sample of Contec HydroPure please contact Contec at infoeu@contecinc.com or by calling +33 (0) 2 97 43 76 90.

*Tested to EN13704 / EN13697

Use biocides safely. Always read the label and product information before use.