

PRODUCT SUMMARY LIST OCT 2016

Product Code	Product Name	Description	Pack Size	Shelf Life	Outer Bags	Endo Tested	Manufacturing Zone
SBT0570IW	Contec Sterile IPA	Sterile 70% Isopropanol in water for injection	8 x 0.5L	3 years	3	Y	Grade A air in Grade C
SBT170IW	Contec Sterile IPA	Sterile 70% Isopropanol in water for injection	6 x 1L	3 years	3	Y	Grade A air in Grade C
SBC570I	Contec Sterile IPA	Sterile 70% Isopropanol in purified water	2 x 5L	3 years	2	Y	Grade A air in Grade C
FBT170I	Contec IPA	Filtered 70% Isopropanol in purified water	6 x 1L	3 years	2	N	Grade A air in Grade C
FBC570I	Contec IPA	Filtered 70% Isopropanol in purified water	2 x 5L	3 years	2	N	Grade A air in Grade C
SBT0570DEW	Contec Sterile Denatured Ethanol	Sterile 70% Denatured Ethanol in water for injection	8 x 0.5L	3 years	3	Y	Grade A air in Grade C
SBT170DEW	Contec Sterile Denatured Ethanol	Sterile 70% Denatured Ethanol in water for injection	6 x 1L	3 years	3	Y	Grade A air in Grade C
SBC570DE	Contec Sterile Denatured Ethanol	Sterile 70% Denatured Ethanol in purified water	2 x 5L	3 years	2	Y	Grade A air in Grade C
FBT170DE	Contec Denatured Ethanol	Filtered 70% Denatured Ethanol in purified water	6 x 1L	3 years	2	N	Grade A air in Grade C
FBC570DE	Contec Denatured Ethanol	Filtered 70% Denatured Ethanol in purified water	2 x 5L	3 years	2	N	Grade A air in Grade C
SBT16HP	Contec Sterile HydroPure	Sterile 6% Stabilised Hydrogen Peroxide in water for injection	6 x 1L	2 years	3	Y	UDAF Grade A in Grade C
SBC56HP	Contec Sterile HydroPure	Sterile 6% Stabilised Hydrogen Peroxide in purified water	2 x 5L	2 years	2	Y	UDAF Grade A in Grade C
FBT16HP	Contec HydroPure	Filtered 6% Stabilised Hydrogen Peroxide in purified water	6 x 1L	1 year	2	N	UDAF Grade A in Grade C
FBC56HP	Contec HydroPure	Filtered 6% Stabilised Hydrogen Peroxide in purified water	2 x 5L	1 year	2	N	UDAF Grade A in Grade C
SBT102PC	Contec Sterile ProChlor	Sterile Stabilised Hypochlorous Acid in purified water	6 x 1L	18 Months	2	N	UDAF Grade A in Grade B
SBC502PC	Contec Sterile ProChlor	Sterile Stabilised Hypochlorous Acid in purified water	2 x 5L	18 Months	2	N	UDAF Grade A in Grade B
FBT102PC	Contec ProChlor	Filtered Stabilised Hypochlorous Acid in purified water	6 x 1L	18 Months	2	N	UDAF Grade A in Grade B
FBC502PC	Contec ProChlor	Filtered Stabilised Hypochlorous Acid in purified water	2 x 5L	18 Months	2	N	UDAF Grade A in Grade B
SBT103NK	Contec Sterile NeutraKlean	Sterile Neutral Detergent in purified water	6 x 1L	2 years	2	N	Grade A Air in Grade C
SBC503NK	Contec Sterile NeutraKlean	Sterile Neutral Detergent in purified water	2 x 5L	2 years	2	N	Grade A Air in Grade C
FBT103NK	Contec NeutraKlean	Filtered Neutral Detergent in purified water (EP)	6 x 1L	2 years	2	N	Grade A Air in Grade C
FBC503NK	Contec NeutraKlean	Filtered Neutral Detergent in purified water (EP)	2 x 5L	2 years	2	N	Grade A Air in Grade C
SSWC102PC	Contec SATWipe ProChlor	Sterile Amplitude Theta with Sterile ProChlor	10 x 100 wipes	18 months	2	N	UDAF Grade A in Grade B
FSWC102PC	Contec SATWipe ProChlor	Amplitude Theta with ProChlor	10 x 100 wipes	18 months	2	N	UDAF Grade A in Grade B
SSWP02PC	Contec SATWipe ProChlor	Sterile Melblown Polypropylene with Sterile ProChlor	10 x 120 wipes	18 months	2	N	UDAF Grade A in Grade B
FSWP02PC	Contec SATWipe ProChlor	Melblown Polypropylene with ProChlor	10 x 120 wipes	18 months	2	N	UDAF Grade A in Grade B